Maxwell

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ELECTRONIC PERSONAL SCALE

This scale uses bioelectric body impedance analysis «BIA» that allows to precisely calculate the percentage of fat, water, muscle and bone tissue. The calculations are based on personal data (height, age, sex) input into the unit memory.

DESCRIPTION

1. Platform
2. Data selection and confirmation sensor button «SET»
3. Display
4. Data selection sensor buttons «◀/▶»
5. Weight measurement units selection button «kg/lb/st»
6. Battery compartment lid
7. Contact pads

RECOMMENDATIONS FOR CUSTOMER ATTENTION!

Before using the unit, read these instructions carefully and keep them for further reference.

Use the unit according to its intended purpose only, as it is stated in this user manual. Mishandling the unit can lead to its breakage and cause harm to the user or damage to his/her property.

Using the excess weight analysis function is TOTALLY contraindicated to pregnant women.

Using the excess weight analysis function is TOTALLY contraindicated to people using pacemakers and other medical devices with built-in batteries.

- Handle your scale with care, as it is a precision instrument, do not subject it to high or low temperatures, high humidity, do not expose it to direct sunlight and do not drop it.
- Use the scale away from heating units.
- Avoid getting of liquids into the scale body, this appliance is not waterproof. Do not leave or use the scale in a room with high humidity (above 80%), make sure that the scale does not come into contact with water or other liquids as this may cause failure of scale indication or its damage.
- If there is a protective film on the weighing platform, remove it to avoid static effect.
- Place the scale on a flat, dry and non-slip surface.
- Do not use the scale on carpets or rugs.
- Do not step on the scale with wet feet; do not use the scale if the weighing platform is damp. You can slip and get injured.
- Be careful while stepping on the scale: step on the scale first with one foot; make sure that the scale is steady and only then step on it with both feet.
- Step on the scale and dispose your weight evenly. Stand still while weighing. Do not jump on the scale surface.
- If you are not planning to use the scale for a long time, remove the battery from the battery compartment.

- Strictly follow the polarity according to the markings when installing the battery.
- Replace the batteries in time.
- This unit is not intended for usage by children under 8 years of age.
- This unit is not intended for usage by people (including children over 8) with physical, neural, mental deviations or with insufficient experience or knowledge. Such persons can use this unit only if they are under supervision of a person who is responsible for their safety and if they are given all the necessary and understandable instructions concerning the safe usage of the unit and information about danger that can be caused by its improper usage.
- Do not leave children unattended not to let them use the unit as a toy.
- For children safety reasons do not leave polyethylene bags, used as a packaging, unattended.

Attention! Do not allow children to play with polyethylene bags or packaging film. **Danger of suffocation!**

- Refer to authorized service centers for maintenance.
- Transport the unit in the original package only.
- Keep the scale away from children in a dry cool place.

THE UNIT IS INTENDED FOR HOUSEHOLD USAGE ONLY

IMPORTANT

Electromagnetic compatibility

The electronic scale can be sensitive to electromagnetic emission of other units located in close proximity (such as mobile phones, portable radio transmitters, radio controllers and microwave ovens). In case of such emission (false or inconsistent data indication on the display), use the scale away from the source of interferences.

USING THE SCALE

After unit transportation or storage at low temperature, it is necessary to keep it for at least two hours at room temperature before switching on.

- Unpack the unit; remove any stickers that can prevent unit operation.
- Clean the weighing platform (1) and the scale body with a soft slightly damp cloth, and then wipe dry.

Battery installation

- Remove the battery compartment lid (6) and insert a “CR2032” battery (supplied with the unit), following the polarity.
- Install the battery compartment lid (6) back to its place.

Note: If there is isolation insert in the battery compartment, open the battery compartment lid (6), remove the isolation insert and install the battery compartment lid (6) back to its place.

ENGLISH

Battery replacement

- The symbol «Lo» indicates that the battery is low, in this case replace the battery.
- Open the battery compartment lid (6), replace the battery with a new CR2032 battery, strictly following the polarity, and close the battery compartment lid (6).
- If you are not going to use the scale for a long time, remove the battery from the battery compartment.

WHAT IS EXCESS WEIGHT?

It is known that excessive weight can negatively affect health. Excessive weight can cause obesity, heart problems, atherosclerosis and osteoporosis. Short weight may also lead to osteoporosis. (Osteoporosis is a systemic skeleton disease characterized by reduction of bone mass per unit volume and impairment of bone tissue structure and leading to increased bone fragility and risk of fracture by minimal trauma or even without it).

It is known that in case of normal weight muscles and body tissues contain more water than fat, and that electric current passes through water better than through fat. Sending weak current through the body, the scale quickly and easily calculates the percentage of fat and water in a human body.

Physical exercises and a moderate diet help to lose excess weight and this scale helps to control your weight.

RECOMMENDATIONS

- Weigh yourself at the same time of the day.
- It is recommended to measure the fat and water percentage in the evening before supper.
- At least 2-3 hours should pass after the meal.
- At least 10-12 hours should pass after intensive physical exercises.
- Percentage of fat in the body depends on different factors and can change during the day.
- When calculating percentage of fat and water of children under 6 years of age, the results may be inaccurate.
- Professional sportsmen, people inclined to edemas and people with nonstandard body proportions can show inaccurate results of calculation.

Note: *Percentage of fat, bone, muscles, and water can be calculated only if personal data (sex, age, height, physical training) are entered correctly and if you step on the scale barefoot.*

SWITCHING ON

- To switch the scale mode on, slightly press on the platform (1), weight measurement units will appear on the display (3).
- To select weight measurement units, press the «kg/lb/st» button (5) consequently. Measurement units are shown on the display (3):
- «kg» - weight in kilograms;
- «lb» - weight in pounds;
- «st» - weight in stones.

- Place the scale on a flat solid surface; if the floor is uneven or covered with a carpet or rug, the weighing accuracy decreases.

WEIGHT MEASUREMENT

- Step on the weighing platform (1).
- The symbols «0.0 kg», «0.0 lb» or «0:0 st» will appear on the display, depending on the selected measurement units.
- Keep your feet parallel to each other and distribute your weight evenly on the platform (1).
- Stay still while weighing. The weight indications on the display (3) will flash during weighing, when the weight is fixed, the indications will light constantly.

Attention!

- **If the symbol «Err» appears on the display (3), it means that the scale is overloaded. Step off the scale immediately to avoid its damage.**
- **The scale will be switched off automatically approximately 20 seconds after weight measurement.**

ENTERING PERSONAL DATA

- Touch the «SET» button (2) twice, one of the user numbers «1», «2».., «12» will be flashing on the display (3).
- Select the user number touching the buttons (4) «◀/▶». You can enter data for 12 users marked by symbols «1», «2».., «12».
- Confirm the user selection touching the «SET» button (2).
- Then the flashing symbol of the user's sex will be shown on display (3). Touching the buttons (4) «◀/▶», select the user's sex and physical training level:

- man;

- man with high physical training level;

- woman;

- woman with high physical training level.

- Confirm your choice by the «SET» button (2).
- Values of age will be flashing on the display (3).
- Set the user's age touching the buttons (4) «◀/▶», the age range is from 6 to 100 years with the step of 1 year. To save the user's age, touch the «SET» button (2).
- Then the height values will be flashing on the display (3). The height range is from 80 to 220 cm with the step of 1 cm. Set the height by touching the buttons (4) «◀/▶». Confirm the set value by touching the «SET» button (2).
- Then the weight measurement units will be flashing on display (3). Touching the buttons (4) «◀/▶», select the necessary weight measurement units (kg/lb/st) and save them by touching the button (2) «SET».

- If it is not necessary to change the weight measurement units, just touch the «SET» button (2).

Note: If during the data entering no actions are taken, the scale will be switched off automatically in approximately 20 seconds.

CALCULATION OF FAT, BONE AND MUSCLE TISSUE AND WATER PERCENTAGE IN A HUMAN BODY

- Touch the «SET» button (2), the symbols «0.0 kg», «0.0 lb» or «0:0 st» will appear on the display (3), depending on the selected measurement units.
- Select your user number by touching the buttons (4) «◀/▶».
- Confirm the user selection touching the «SET» button (2).
- Touching the «SET» button (2) repeatedly, check the user's previously saved personal data (sex, physical training, age, height).
- Then you'll see a weight measurement unit symbol flashing on the display (3), confirm the choice of the weight measurement unit by touching the «SET» button (2).
- Step with your bare feet on the scale so that your feet touch the metal contacts on the platform (1). Stay still during the weight measurement and wait till the weight indications appear.
- Keep standing on the scale. Symbols «- - -» on the display (3) indicate that the percentage of fat, water, bone and muscles tissue is being calculated.
- After calculation the following data will be shown on the display (3):
 - the user's weight in selected measurement units (kg/lb/st);

FAT

- fat percentage,

WAT

- water percentage,

MUS

- muscle tissue percentage,

BON

- bone tissue percentage,

KCAL – recommended level of daily calorie intake,

«-» - low weight,

«- - -» - normal weight,

«- - - - -» - excessive weight,

«- - - - - - -» - obesity.

- Calculation results will be shown on the display (3) three times and then the scale will be switched off automatically.

Attention!

- If the symbol «Err» appears on the display (3), it means that the scale is overloaded. Step off the scale immediately to avoid its damage.
- The symbol «Err-2» shows that an error occurred during calculation. Step off the scale and repeat weighing.

Note: The data concerning the percentage of fat and water may be obtained only if you have entered the user data and stand on the scale barefoot.

Change of fat percentage in a human body during the day

The indications of fat percentage in a human body are affected by such factors as lifestyle, physical exercises, health condition and amount of liquid in the body.

You should take into account the following factors:

- the fat percentage indications may be less if measured after taking a bath or a shower,
- the fat percentage indications may be higher if measured after having a meal or drinking,
- percentage of liquid in a human body may decrease during and after illness and after physical exercises and thus change the fat percentage indications.

Water percentage in a human body

Water percentage in a human body is the percentage of liquid in the human body to total weight. Water has primary importance for human tissues and organs, since 50-65% of our body consist of water. With this scale you can control the percentage of water in your body and correct this important characteristic.

What is muscle bulk?

Muscle bulk value shows the weight of locomotor system muscles, plain muscles (heart and digestive system) and water in these muscles.

Control of muscle bulk is especially important for people with excessive weight and obesity. Weight loss may lead to decrease of fat tissues, but only if you gain muscle bulk at the same time.

What is bone mass?

Bone tissues are the base of human skeleton; they are solid and resistant to external actions, protecting internal organs, and together with the muscle tissue providing the ability of a human to move.

Bone mass measurement is based on calculating the minerals mass in a human body (calcium or other matters). The research showed that development of muscle tissues makes the skeleton stronger. This fact should be taken into account when making up a diet and calculating physical activity for strong and healthy bones. Bone mass practically doesn't change; however, some decrease is possible under many factors (weight, height, age, sex). That's why this manual does not contain standard values and recommendations concerning the bone tissue percentage.

Recommended level of daily calorie intake

The level of daily calorie intake is the amount of energy required by the human body at absolute rest for maintaining its vital functions (for instance, if you stay in bed for 24 hours). The required amount of calories significantly depends on the weight, height and age of a person. This amount of calories is required by the human body and it should be regained by means of food. If you receive less calories for a long time, it may harm your health.

Normal percentage of water in the body is 45-60% for women and 50-65% for men.

ENGLISH

Man	Age	% of fat			% of water			% of muscles		
		Low	Medium	High	Low	Medium	High	Low	Medium	High
	20-39	5.0-7.9	8.0-19.0	19.1-45.0	38.0-54.9	55.0-70.0	70.1-85.0	45.0-72.9	73.0-82.0	82.1-95.0
	40-59	5.0-10.9	11.0-22.0	22.1-45.0	38.0-49.9	50.0-68.0	68.1-85.0	45.0-69.9	70.0-79.0	79.1-95.0
	60-79	5.0-12.9	13.0-25.0	25.1-45.0	38.0-49.9	50.0-65.0	65.1-85.0	45.0-66.9	67.0-77.0	77.1-95.0
Woman		% of fat			% of water			% of muscles		
	20-39	5.0-20.9	21.0-33.0	33.1-45.0	38.0-47.9	48.0-65.0	65.1-85.0	45.0-67.9	68.0-77.0	77.1-95.0
	40-59	5.0-22.9	23.0-35.0	35.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0	45.0-64.9	65.0-75.0	75.1-95.0
	60-79	5.0-23.9	24.0-36.0	36.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0	45.0-61.9	62.0-72.0	72.1-95.0
Man taking exercises up to 10 hours a week		% of fat			% of water			% of muscles		
	20-39	5.0-5.9	6.0-12.0	12.1-45.0	38.0-54.9	55.0-71.0	71.1-85.0	45.0-74.9	75.0-85.0	85.1-95.0
	40-59	5.0-5.9	6.0-13.0	13.1-45.0	38.0-49.9	50.0-69.0	69.1-85.0	45.0-72.9	73.0-83.0	83.1-95.0
	60-79	5.0-7.9	8.0-13.0	13.1-45.0	38.0-49.9	50.0-67.0	67.1-85.0	45.0-69.9	70.0-80.0	80.1-95.0
Woman taking exercises up to 10 hours a week		% of fat			% of water			% of muscles		
	20-39	5.0-13.9	14.0-18.0	18.1-45.0	38.0-47.9	48.0-66.0	66.1-85.0	45.0-67.9	68.0-80.0	80.1-95.0
	40-59	5.0-15.9	16.0-20.0	20.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0	45.0-65.9	66.0-78.0	78.1-95.0
	60-79	5.0-16.9	17.0-20.0	20.1-45.0	38.0-44.9	45.0-65.0	65.1-85.0	45.0-64.9	65.0-76.0	76.1-95.0

Note: The data in this table are given for reference.

Analysis of weighing results

As a rule, short-term changes of weight within several days are the result of minor loss of liquid in your body tissues. Analyze the results of regular weight measurements for a long time. Evaluation of the results is based on changes of total weight, fat and water percentage in body tissues and depends on the duration of the period when the measurements are made.

Short-term indication changes during a day can be significantly different from the indication changes within a week or a month.

As a rule, short-term indication changes are connected with percentage of liquid in body tissues, and medium and long-term changes are caused by changes in the amount of fat and muscle bulk.

- Short-term weight loss and simultaneous growth or constant amount of fat in body tissues may be caused by decreasing of liquid amount as a result of physical exercises, visiting a sauna or at quick weight loss.
- Medium-term increase of weight, decrease of fat amount or constant fat amount in body tissues may be caused by increasing of muscle bulk.
- If you notice simultaneous decrease of total weight and fat percentage, it means that your diet works, i.e. you lose the fat content in your body tissues.
- The best results are achieved if you combine diet with physical exercises. In this case you can observe the medium-term increase of muscle mass in body tissues.

CLEANING AND CARE

1. Clean the scale with a soft, slightly damp cloth and then wipe it dry.

2. Never use abrasives or solvents.
3. Never immerse the scale into water or other liquids.
4. Remove the battery from the battery compartment before taking the unit away for storage.
5. Keep the scale away from children in a dry cool place.

DELIVERY SET

- Scale – 1 pc.
- CR2032 battery – 1 pc.
- Instruction manual – 1 pc.

TECHNICAL SPECIFICATIONS

Power supply: one CR2032 battery, 3 V
Maximal weight: 180 kg

The manufacturer preserves the right to change the design and the specifications of the unit without a preliminary notification.

Unit operating life is 3 years

Guarantee

Details regarding guarantee conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

 This product conforms to the EMC-Requirements as laid down by the Council Directive 2004/108/EC and to the Low Voltage Regulation (2006/95/EC)

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© ООО ГОЛДЕР-ЭЛЕКТРОНИКС, 2013
© GOLDER-ELECTRONICS LLC, 2013